

IDDSI adoption and implementation is happening across Canada. In order to help facilitate this initiative, we have mapped many of our existing products to determine their IDDSI label and level using IDDSI approved audit tools.

PRODUCT	 PUREED	 MINCED & MOIST	 SOFT & BITE-SIZED	 EASY TO CHEW
Basa fillet, 4oz, IQF (SF03102)				
Halibut fillet, 4oz, IQF (SF03504)				
Pollock fillet, 4-6oz, IQF (SF00502, SF01501)				
Salmon fillet, 6oz & 8oz, IQF (SF04641, SF0464)				
Tilapia fillet, 3-5oz, IQF (SF03138)				
PouchPak™ Tuna (SF00965, SF00990)				
Crab Cakes, 3oz, IQF (SF04151)	 *			

*DID NOT MEET IDDSI CRITERIA FOR LEVEL 4 PUREED.

IDDSI Mapping Summary

IDDSI Level	Items Meeting Criteria	Preparation	Particle Size	Critical Tests Used	Meets Criteria At		
					Time of Service	15 min After Serving	30 min After Serving
4 PUREED	Basa, Halibut, Pollock, Salmon, Tilapia, Tuna	Samples were prepared, mechanically processed with vegetable broth and commercial thickener and tested	N/A	Appearance, Fork Drip Test, Spoon Tilt Test	✓	✓	✓
5 MINCED & MOIST	Basa, Crab Cake, Halibut, Pollock, Salmon, Tilapia, Tuna	Samples were prepared, mashed with fork and tested	Food intended to mimic “chewed bolus”, equal/less than 4mm x 15mm (adults)	Appearance, Fork Drip Test, Spoon Tilt Test OR Finger Test if these are not available	✓	N/A	N/A
6 SOFT & BITE-SIZED	Basa, Crab Cake, Halibut, Pollock, Salmon, Tilapia, Tuna	Samples were prepared, cut into smaller pieces and tested	Food intended to mimic a “bite of food”, equal/less than 15mm x 15mm (adults)	Appearance, Fork Drip Test, Spoon Tilt Test OR Finger Test if these are not available	✓	N/A	N/A
7 EASY TO CHEW	Basa, Crab Cake, Halibut, Pollock, Salmon, Tilapia, Tuna	Samples were prepared, fork cut and tested	Fork cut to 15mm x 15mm pieces	Fork/Spoon Pressure Test OR Finger Test if these are not available	✓	N/A	N/A

DISCLAIMER: PLEASE NOTE THAT RESULTS MAY VARY DEPENDENT ON METHODS OF PREPARATION SPECIFIC TO YOUR OPERATION.